

Attorney Docket No. 0225-0068.30

Patent

**Transmittal of Utility Patent Application for Filing**

*Certification Under 37 C.F.R. §1.10 (if applicable)*

EL 530 372 242 US  
"Express Mail" Label Number

July 11, 2001  
Date of Deposit

I hereby certify that this application, and any other documents referred to as enclosed herein are being deposited in an envelope with the United States Postal Service "Express Mail Post Office to Addressee" service under 37 CFR §1.10 on the date indicated above and addressed to the Assistant Commissioner for Patents, Washington, D.C. 20231

**Jennifer L. Mahoney**

(Print Name of Person Mailing Application)

(Signature of Person Mailing Application)



**TAG CLEAVAGE FOR DETECTION OF NUCLEIC ACIDS**

This application claims the benefit of U.S. Provisional Application No. 60/217,624 filed July 11, 2000, which is incorporated herewith by reference in its entirety.

**Technical Field**

The field of this invention is the detection of nucleic acids, particularly in a multiplexed protocol.

**Background**

As the elucidation of the human genome and other genomes nears, there is increasing interest in how the information may be used to enhance medical practice. There is increasing evidence that the phenotype of an individual plays a crucial role in the health and well being of an individual. The ability to ward off disease, respond to infection, be resistant to cancer, neurological diseases, autoimmune diseases, and the like, appear to be intimately related to the phenotype. Inherited diseases, associated with mutations in a gene are well known and frequently occur at various phases during an individual's lifetime. There is also interest in relating phenotype to a repertoire of individual variations in the genotype. This requires that one have knowledge of the frequency with which a nucleotide is present at a specific site in a chromosome.

15

20

In many situations, particularly with single nucleotide polymorphisms ("snps"), there is an interest in making a plurality of determinations simultaneously. To this end, one wishes to have techniques that do not result in interference between the different components of the sample, provide for amplification of the nucleic acid sequences that are of interest without significant amplification of other nucleic acid sequences, permit miniaturization to minimize the amount of sample and reagents required for the determination, and the like.

There are a number of techniques for determining nucleic acid sequences and snps. One technique employs two nucleic acid sequences, a primer and a labeled probe, where the primer is extended in the presence of a cleavase. Extension of the primer results in cleaving the label from the probe. The Taqman® process employs a label involving a spaced apart fluorescer and quencher. When the correct sequence is present, the link between the fluorescer and quencher is cleaved, resulting in an increase in fluorescence. In many cases, the methods require the use of the polymerase chain reaction ("PCR"). While PCR is efficient in amplifying a few sequences in the same reaction, once the number of sequences to be amplified exceeds a relatively low level, the efficiency is substantially reduced.

There is, therefore, substantial interest in having protocols that are robust, allow for amplification of a large number of members of a sample without interference, and provide signals that are readily detected and can be individually identified.

#### Brief Description of the Prior Art

References of interest include WO97/740-28 AND 97/00967.

25

#### SUMMARY OF THE INVENTION

Methods and compositions are provided for detecting nucleic acids, either a sequence of at least two nucleotides or a single nucleotide. The method employs a primer and a tagged probe, where the primer is bonded to a member of a bond cleaving system and the probe comprises a detectable tag bonded to a link susceptible to cleavage by said bond cleaving system. The primer and probe bind to a target nucleic acid, whereby the detectable tag is released. Desirably, the hybridization is performed

at a temperature near the melting temperature of the complex of the probe with the target nucleic acid, so that the binding of the probe can be cycled with a plurality of detectable tags released for a single target molecule. The detectable tags may then be analyzed to determine the presence of the target nucleic acid. By employing detectable tags that can be individually analyzed, a plurality of nucleic acid moieties may be detected in the same sample medium.

#### DESCRIPTION OF THE SPECIFIC EMBODIMENTS

Nucleic acid analysis is performed with a combination of a primer and a probe.

10 The primer is characterized by having a nucleic acid binding sequence that binds to a target nucleic acid. Bonded to the sequence is a member of a bond-cleaving system. The probe comprises a nucleic acid binding sequence that hybridizes to a sequence of the target nucleic acid in proximity to the sequence to which the primer hybridizes. Bonded to the probe sequence is a detectable tag, which is bonded to a link susceptible to cleavage by the bond-cleaving system. Desirably, the melting temperature of the complex formed by the probe with the target nucleic acid is lower than the melting temperature of the complex formed by the primer with the target nucleic acid and the process is performed at a temperature at which the primer hybridization is substantially stable, while the primer hybridization, particularly after cleavage by the bond-cleaving system is substantially labile.

20 The method comprises combining under hybridizing conditions, the sample containing the target nucleic acid, and one or more pairs of primers and probes in a liquid medium. The temperature will depend upon whether cycling of the probe is intended, in which case the temperature will be between the melting temperature of the complexes of the template and probe, as influenced by the hybridization between the primer and probe, usually not more than about 10°C greater than the melting temperature of the complex with the probe, conveniently not more than about 5°C greater, and below the melting temperature of the primer and template. After sufficient time for hybridization and cleavage to occur, the liberated tags are analyzed in accordance with the nature of the tag.

The primer has two primary components joined by a link. The first component is a sequence that can hybridize to a target nucleic acid sequence. This sequence will have at least about 12 units (where a unit intends a base), more usually at least about 18 units and not more than about 100 units, usually not more than about 60 units and 5 preferably not more than about 36 units. There will be a sufficient number of units to provide the desired affinity and specificity for the target and difference in binding affinity between the primer and probe. For the most part, the sequence will comprise naturally occurring bases, pyrimidines and purines, linked by natural or unnatural linkages, e.g. phosphate ribose or deoxyribose esters, phosphate ribose or deoxyribose thioesters, 10 phosphate ribose or deoxyribose amides, amino acids, particularly glycine or alanine amides, where the base may be pendent from the amino nitrogen, or the like.

The primer also has an effector moiety, which is normally present as an arm that is not bound to the target sequence and comprises a member of a susceptible bond cleaving system. The effector moiety will usually be capable of cleaving a linkage 15 between the detectable tag and the probe sequence in conjunction with other agents or may provide a complex that results in providing a labile linkage in the probe. The probe sequence may be preprepared, so as to have the detectable tag bonded to the nucleic acid sequence when added to the reaction mixture or result from the addition of one or 20 more nucleotides, usually one nucleotide, where the added nucleotide(s) includes the detectable tag linked to the nucleotide by means of a cleavable linker. The probe will be characterized by having a sequence, which hybridizes to a sequence of target nucleic acid in proximity to the sequence to which the primer binds. Usually, the probe and 25 primer will be separated by fewer than six nucleotides along the chain of the target sequence, usually not more than three nucleotides and may be contiguous (0 – 6; 0 – 3). As indicated, the probe will ultimately have proximal to its terminus, either 3' or 5', the detectable tag linked through a cleavable linkage to the target nucleic acid binding sequence.

The primer will usually have at least the same number of nucleotides binding to 30 the target sequence as the probe and usually more, generally at least about 10% more, more usually at least about 50% more and not more than about 5x, usually not more than about 3x the number of nucleotides. Generally the primer will have at least about

15 bases, more usually at least about 18 bases and not more than about 75 bases, generally not more than about 36 bases, binding to the target sequence. While the total number of bases above the indicated number for the primer is not critical, for the most part there will be no advantage to having a greater number of bases, and the larger the 5 number of bases, the greater the cost.

The probe will also have two regions, the region binding to the target sequence and an arm that is not bound to the target sequence and comprises the cleavable linkage. The probe will have at least about 12 bases, usually at least about 15 bases and usually not more than about 36 bases binding to the target sequence. Desirably 10 the probe will add to the specificity of the assay, when taken together with the primer, but will also allow for cycling. By performing the assay at a temperature between the melting temperature of the complex between the probe and the target nucleic acid and the melting temperature of the complex between the primer and the target nucleic acid, the primer will be substantially retained bound to the target nucleic acid, while the probe will come on and off. During the probe's residence time bound to the target nucleic acid, the opportunity exists for the cleavable link to be cleaved by the effector while bound to the target nucleic acid. By providing for some binding energy between the unbound arms of the primer and probe, upon cleavage of the labile linkage, the melting temperature of the probe will be lowered allowing for greater release. In this manner, 15 one may amplify the number of detectable tags that are released in relation to the number of target sequences present in the sample.

The bond-cleaving system may involve a variety of different components. One or more of the components may be on the arm of the primer. The primer associated component(s) may include an effector which may be an active participant in being 25 involved in being a chemical component of a chemical reaction or may be a passive participant, providing an environment recognized by an enzyme, where the enzyme will act as a hydrolase, or both. The effector normally will cooperate with another agent that may be a chemical or physical agent, such as electromagnetic radiation or the like. Therefore, the effector will normally be unchanged or cycled. The effector may be 30 organic or inorganic, if inorganic, usually metalloorganic, will usually require additional agents, either physical or chemical. Included among the effectors are enzymes,

15  
20  
25  
30  
35  
40  
45  
50  
55  
60  
65  
70  
75  
80  
85  
90  
95  
100  
105  
110  
115  
120  
125  
130  
135  
140  
145  
150  
155  
160  
165  
170  
175  
180  
185  
190  
195  
200  
205  
210  
215  
220  
225  
230  
235  
240  
245  
250  
255  
260  
265  
270  
275  
280  
285  
290  
295  
300  
305  
310  
315  
320  
325  
330  
335  
340  
345  
350  
355  
360  
365  
370  
375  
380  
385  
390  
395  
400  
405  
410  
415  
420  
425  
430  
435  
440  
445  
450  
455  
460  
465  
470  
475  
480  
485  
490  
495  
500  
505  
510  
515  
520  
525  
530  
535  
540  
545  
550  
555  
560  
565  
570  
575  
580  
585  
590  
595  
600  
605  
610  
615  
620  
625  
630  
635  
640  
645  
650  
655  
660  
665  
670  
675  
680  
685  
690  
695  
700  
705  
710  
715  
720  
725  
730  
735  
740  
745  
750  
755  
760  
765  
770  
775  
780  
785  
790  
795  
800  
805  
810  
815  
820  
825  
830  
835  
840  
845  
850  
855  
860  
865  
870  
875  
880  
885  
890  
895  
900  
905  
910  
915  
920  
925  
930  
935  
940  
945  
950  
955  
960  
965  
970  
975  
980  
985  
990  
995  
1000  
1005  
1010  
1015  
1020  
1025  
1030  
1035  
1040  
1045  
1050  
1055  
1060  
1065  
1070  
1075  
1080  
1085  
1090  
1095  
1100  
1105  
1110  
1115  
1120  
1125  
1130  
1135  
1140  
1145  
1150  
1155  
1160  
1165  
1170  
1175  
1180  
1185  
1190  
1195  
1200  
1205  
1210  
1215  
1220  
1225  
1230  
1235  
1240  
1245  
1250  
1255  
1260  
1265  
1270  
1275  
1280  
1285  
1290  
1295  
1300  
1305  
1310  
1315  
1320  
1325  
1330  
1335  
1340  
1345  
1350  
1355  
1360  
1365  
1370  
1375  
1380  
1385  
1390  
1395  
1400  
1405  
1410  
1415  
1420  
1425  
1430  
1435  
1440  
1445  
1450  
1455  
1460  
1465  
1470  
1475  
1480  
1485  
1490  
1495  
1500  
1505  
1510  
1515  
1520  
1525  
1530  
1535  
1540  
1545  
1550  
1555  
1560  
1565  
1570  
1575  
1580  
1585  
1590  
1595  
1600  
1605  
1610  
1615  
1620  
1625  
1630  
1635  
1640  
1645  
1650  
1655  
1660  
1665  
1670  
1675  
1680  
1685  
1690  
1695  
1700  
1705  
1710  
1715  
1720  
1725  
1730  
1735  
1740  
1745  
1750  
1755  
1760  
1765  
1770  
1775  
1780  
1785  
1790  
1795  
1800  
1805  
1810  
1815  
1820  
1825  
1830  
1835  
1840  
1845  
1850  
1855  
1860  
1865  
1870  
1875  
1880  
1885  
1890  
1895  
1900  
1905  
1910  
1915  
1920  
1925  
1930  
1935  
1940  
1945  
1950  
1955  
1960  
1965  
1970  
1975  
1980  
1985  
1990  
1995  
2000  
2005  
2010  
2015  
2020  
2025  
2030  
2035  
2040  
2045  
2050  
2055  
2060  
2065  
2070  
2075  
2080  
2085  
2090  
2095  
2100  
2105  
2110  
2115  
2120  
2125  
2130  
2135  
2140  
2145  
2150  
2155  
2160  
2165  
2170  
2175  
2180  
2185  
2190  
2195  
2200  
2205  
2210  
2215  
2220  
2225  
2230  
2235  
2240  
2245  
2250  
2255  
2260  
2265  
2270  
2275  
2280  
2285  
2290  
2295  
2300  
2305  
2310  
2315  
2320  
2325  
2330  
2335  
2340  
2345  
2350  
2355  
2360  
2365  
2370  
2375  
2380  
2385  
2390  
2395  
2400  
2405  
2410  
2415  
2420  
2425  
2430  
2435  
2440  
2445  
2450  
2455  
2460  
2465  
2470  
2475  
2480  
2485  
2490  
2495  
2500  
2505  
2510  
2515  
2520  
2525  
2530  
2535  
2540  
2545  
2550  
2555  
2560  
2565  
2570  
2575  
2580  
2585  
2590  
2595  
2600  
2605  
2610  
2615  
2620  
2625  
2630  
2635  
2640  
2645  
2650  
2655  
2660  
2665  
2670  
2675  
2680  
2685  
2690  
2695  
2700  
2705  
2710  
2715  
2720  
2725  
2730  
2735  
2740  
2745  
2750  
2755  
2760  
2765  
2770  
2775  
2780  
2785  
2790  
2795  
2800  
2805  
2810  
2815  
2820  
2825  
2830  
2835  
2840  
2845  
2850  
2855  
2860  
2865  
2870  
2875  
2880  
2885  
2890  
2895  
2900  
2905  
2910  
2915  
2920  
2925  
2930  
2935  
2940  
2945  
2950  
2955  
2960  
2965  
2970  
2975  
2980  
2985  
2990  
2995  
3000  
3005  
3010  
3015  
3020  
3025  
3030  
3035  
3040  
3045  
3050  
3055  
3060  
3065  
3070  
3075  
3080  
3085  
3090  
3095  
3100  
3105  
3110  
3115  
3120  
3125  
3130  
3135  
3140  
3145  
3150  
3155  
3160  
3165  
3170  
3175  
3180  
3185  
3190  
3195  
3200  
3205  
3210  
3215  
3220  
3225  
3230  
3235  
3240  
3245  
3250  
3255  
3260  
3265  
3270  
3275  
3280  
3285  
3290  
3295  
3300  
3305  
3310  
3315  
3320  
3325  
3330  
3335  
3340  
3345  
3350  
3355  
3360  
3365  
3370  
3375  
3380  
3385  
3390  
3395  
3400  
3405  
3410  
3415  
3420  
3425  
3430  
3435  
3440  
3445  
3450  
3455  
3460  
3465  
3470  
3475  
3480  
3485  
3490  
3495  
3500  
3505  
3510  
3515  
3520  
3525  
3530  
3535  
3540  
3545  
3550  
3555  
3560  
3565  
3570  
3575  
3580  
3585  
3590  
3595  
3600  
3605  
3610  
3615  
3620  
3625  
3630  
3635  
3640  
3645  
3650  
3655  
3660  
3665  
3670  
3675  
3680  
3685  
3690  
3695  
3700  
3705  
3710  
3715  
3720  
3725  
3730  
3735  
3740  
3745  
3750  
3755  
3760  
3765  
3770  
3775  
3780  
3785  
3790  
3795  
3800  
3805  
3810  
3815  
3820  
3825  
3830  
3835  
3840  
3845  
3850  
3855  
3860  
3865  
3870  
3875  
3880  
3885  
3890  
3895  
3900  
3905  
3910  
3915  
3920  
3925  
3930  
3935  
3940  
3945  
3950  
3955  
3960  
3965  
3970  
3975  
3980  
3985  
3990  
3995  
4000  
4005  
4010  
4015  
4020  
4025  
4030  
4035  
4040  
4045  
4050  
4055  
4060  
4065  
4070  
4075  
4080  
4085  
4090  
4095  
4100  
4105  
4110  
4115  
4120  
4125  
4130  
4135  
4140  
4145  
4150  
4155  
4160  
4165  
4170  
4175  
4180  
4185  
4190  
4195  
4200  
4205  
4210  
4215  
4220  
4225  
4230  
4235  
4240  
4245  
4250  
4255  
4260  
4265  
4270  
4275  
4280  
4285  
4290  
4295  
4300  
4305  
4310  
4315  
4320  
4325  
4330  
4335  
4340  
4345  
4350  
4355  
4360  
4365  
4370  
4375  
4380  
4385  
4390  
4395  
4400  
4405  
4410  
4415  
4420  
4425  
4430  
4435  
4440  
4445  
4450  
4455  
4460  
4465  
4470  
4475  
4480  
4485  
4490  
4495  
4500  
4505  
4510  
4515  
4520  
4525  
4530  
4535  
4540  
4545  
4550  
4555  
4560  
4565  
4570  
4575  
4580  
4585  
4590  
4595  
4600  
4605  
4610  
4615  
4620  
4625  
4630  
4635  
4640  
4645  
4650  
4655  
4660  
4665  
4670  
4675  
4680  
4685  
4690  
4695  
4700  
4705  
4710  
4715  
4720  
4725  
4730  
4735  
4740  
4745  
4750  
4755  
4760  
4765  
4770  
4775  
4780  
4785  
4790  
4795  
4800  
4805  
4810  
4815  
4820  
4825  
4830  
4835  
4840  
4845  
4850  
4855  
4860  
4865  
4870  
4875  
4880  
4885  
4890  
4895  
4900  
4905  
4910  
4915  
4920  
4925  
4930  
4935  
4940  
4945  
4950  
4955  
4960  
4965  
4970  
4975  
4980  
4985  
4990  
4995  
5000  
5005  
5010  
5015  
5020  
5025  
5030  
5035  
5040  
5045  
5050  
5055  
5060  
5065  
5070  
5075  
5080  
5085  
5090  
5095  
5100  
5105  
5110  
5115  
5120  
5125  
5130  
5135  
5140  
5145  
5150  
5155  
5160  
5165  
5170  
5175  
5180  
5185  
5190  
5195  
5200  
5205  
5210  
5215  
5220  
5225  
5230  
5235  
5240  
5245  
5250  
5255  
5260  
5265  
5270  
5275  
5280  
5285  
5290  
5295  
5300  
5305  
5310  
5315  
5320  
5325  
5330  
5335  
5340  
5345  
5350  
5355  
5360  
5365  
5370  
5375  
5380  
5385  
5390  
5395  
5400  
5405  
5410  
5415  
5420  
5425  
5430  
5435  
5440  
5445  
5450  
5455  
5460  
5465  
5470  
5475  
5480  
5485  
5490  
5495  
5500  
5505  
5510  
5515  
5520  
5525  
5530  
5535  
5540  
5545  
5550  
5555  
5560  
5565  
5570  
5575  
5580  
5585  
5590  
5595  
5600  
5605  
5610  
5615  
5620  
5625  
5630  
5635  
5640  
5645  
5650  
5655  
5660  
5665  
5670  
5675  
5680  
5685  
5690  
5695  
5700  
5705  
5710  
5715  
5720  
5725  
5730  
5735  
5740  
5745  
5750  
5755  
5760  
5765  
5770  
5775  
5780  
5785  
5790  
5795  
5800  
5805  
5810  
5815  
5820  
5825  
5830  
5835  
5840  
5845  
5850  
5855  
5860  
5865  
5870  
5875  
5880  
5885  
5890  
5895  
5900  
5905  
5910  
5915  
5920  
5925  
5930  
5935  
5940  
5945  
5950  
5955  
5960  
5965  
5970  
5975  
5980  
5985  
5990  
5995  
6000  
6005  
6010  
6015  
6020  
6025  
6030  
6035  
6040  
6045  
6050  
6055  
6060  
6065  
6070  
6075  
6080  
6085  
6090  
6095  
6100  
6105  
6110  
6115  
6120  
6125  
6130  
6135  
6140  
6145  
6150  
6155  
6160  
6165  
6170  
6175  
6180  
6185  
6190  
6195  
6200  
6205  
6210  
6215  
6220  
6225  
6230  
6235  
6240  
6245  
6250  
6255  
6260  
6265  
6270  
6275  
6280  
6285  
6290  
6295  
6300  
6305  
6310  
6315  
6320  
6325  
6330  
6335  
6340  
6345  
6350  
6355  
6360  
6365  
6370  
6375  
6380  
6385  
6390  
6395  
6400  
6405  
6410  
6415  
6420  
6425  
6430  
6435  
6440  
6445  
6450  
6455  
6460  
6465  
6470  
6475  
6480  
6485  
6490  
6495  
6500  
6505  
6510  
6515  
6520  
6525  
6530  
6535  
6540  
6545  
6550  
6555  
6560  
6565  
6570  
6575  
6580  
6585  
6590  
6595  
6600  
6605  
6610  
6615  
6620  
6625  
6630  
6635  
6640  
6645  
6650  
6655  
6660  
6665  
6670  
6675  
6680  
6685  
6690  
6695  
6700  
6705  
6710  
6715  
6720  
6725  
6730  
6735  
6740  
6745  
6750  
6755  
6760  
6765  
6770  
6775  
6780  
6785  
6790  
6795  
6800  
6805  
6810  
6815  
6820  
6825  
6830  
6835  
6840  
6845  
6850  
6855  
6860  
6865  
6870  
6875  
6880  
6885  
6890  
6895  
6900  
6905  
6910  
6915  
6920  
6925  
6930  
6935  
6940  
6945  
6950  
6955  
6960  
6965  
6970  
6975  
6980  
6985  
6990  
6995  
7000  
7005  
7010  
7015  
7020  
7025  
7030  
7035  
7040  
7045  
7050  
7055  
7060  
7065  
7070  
7075  
7080  
7085  
7090  
7095  
7100  
7105  
7110  
7115  
7120  
7125  
7130  
7135  
7140  
7145  
7150  
7155  
7160  
7165  
7170  
7175  
7180  
7185  
7190  
7195  
7200  
7205  
7210  
7215  
7220  
7225  
7230  
7235  
7240  
7245  
7250  
7255  
7260  
7265  
7270  
7275  
7280  
7285  
7290  
7295  
7300  
7305  
7310  
7315  
7320  
7325  
7330  
7335  
7340  
7345  
7350  
7355  
7360  
7365  
7370  
7375  
7380  
7385  
7390  
7395  
7400  
7405  
7410  
7415  
7420  
7425  
7430  
7435  
7440  
7445  
7450  
7455  
7460  
7465  
7470  
7475  
7480  
7485  
7490  
7495  
7500  
7505  
7510  
7515  
7520  
7525  
7530  
7535  
7540  
7545  
7550  
7555  
7560  
7565  
7570  
7575  
7580  
7585  
7590  
7595  
7600  
7605  
7610  
7615  
7620  
7625  
7630  
7635  
7640  
7645  
7650  
7655  
7660  
7665  
7670  
7675  
7680  
7685  
7690  
7695  
7700  
7705  
7710  
7715  
7720  
7725  
7730  
7735  
7740  
7745  
7750  
7755  
7760  
7765  
7770  
7775  
7780  
7785  
7790  
7795  
7800  
7805  
7810  
7815  
7820  
7825  
7830  
7835  
7840  
7845  
7850  
7855  
7860  
7865  
7870  
7875  
7880  
7885  
7890  
7895  
7900  
7905  
7910  
7915  
7920  
7925  
7930  
7935  
7940  
7945  
7950  
7955  
7960  
7965  
7970  
7975  
7980  
7985  
7990  
7995  
8000  
8005  
8010  
8015  
8020  
8025  
8030  
8035  
8040  
8045  
8050  
8055  
8060  
8065  
8070  
8075  
8080  
8085  
8090  
8095  
8100  
8105  
8110  
8115  
8120  
8125  
8130  
8135  
8140  
8145  
8150  
8155  
8160  
8165  
8170  
8175  
8180  
8185  
8190  
8195  
8200  
8205  
8210  
8215  
8220  
8225  
8230  
8235  
8240  
8245  
8250  
8255  
8260  
8265  
8270  
8275  
8280  
8285  
8290  
8295  
8300  
8305  
8310  
8315  
8320  
8325  
8330  
8335  
8340  
8345  
8350  
8355  
8360  
8365  
8370  
8375  
8380  
8385  
8390  
8395  
8400  
8405  
8410  
8415  
8420  
8425  
8430  
8435  
8440  
8445  
8450  
8455  
8460  
8465  
8470  
8475  
8480  
8485  
8490  
8495  
8500  
8505  
8510  
8515  
8520  
8525  
8530  
8535  
8540  
8545  
8550  
8555  
8560  
8565  
8570  
8575  
8580  
8585  
8590  
8595  
8600  
8605  
8610  
8615  
8620  
8625  
8630  
8635  
8640  
8645  
8650  
8655  
8660  
8665  
8670  
8675  
8680  
8685  
8690  
8695  
8700  
8705  
8710  
8715  
8720  
8725  
8730  
8735  
8740  
8745  
8750  
8755  
8760  
8765  
8770  
8775  
8780  
8785  
8790  
8795  
8800  
8805  
8810  
8815  
8820  
8825  
8830  
8835  
8840  
8845  
8850  
8855  
8860  
8865  
8870  
8875  
8880  
8885  
8890  
8895  
8900  
8905  
8910  
8915  
8920  
8925  
8930  
8935  
8940  
8945  
8950  
8955  
8960  
8965  
8970  
8975  
8980  
8985  
8990  
8995  
9000  
9005  
9010  
9015  
9020  
9025  
9030  
9035  
9040  
9045  
9050  
9055  
9060  
9065  
9070  
9075  
9080  
9085  
9090  
9095  
9100  
9

ribozymes, metalloorganic compounds, electromagnetic radiation sensitizers, nucleic acid sequences, or other agents that can be cycled and provide, directly or indirectly, for cleavage of a linkage.

Various enzymes may be used, where the enzymes can cleave one strand of a hybrid DNA-RNA strand, a uracil, a glycol linkage, where the glycol linkage may take other forms than a ribose or deoxyribose diphosphate. The enzymes may be bound to the primer, covalently or non-covalently, or free in solution. A number of enzymes are known that selectively cleave one of the strands of the hybrid complex. These enzymes include RNase H, Rnase H1 and like enzymes from a variety of species, as described in Cerritelli and Crouch, 1998 Genomics 53, 300-307, and references cited therein.

Enzymes that cleave at apurinic and apyrimidinic sites in dsDNA include AP endonucleases, as described in Doetsch and Cunningham, 1990 Mutation Research 236, 173-201; and Levin and Demple, 1990 Nucleic Acids Res. 18, 5069-5075.

Enzymes are available that cleave a DNA strand comprising a uracil, such as uracil-DNA glycosylase, as described in Weiss *et al.*, 1983 Biochemistry 22, 4501-4507. Restriction enzymes that cleave at only one strand include site specific nickases, as described in Abdurashitov, *et al.*, 1996 Molecular Biology 30, 754-758. Other enzymes will cleave at a loop, as a result of a mismatch forming a loop out. Such enzymes include resolvases and specific endonucleases, as described in Tito, jr., *et al.*, 1998 Clin. Chem. 44, 731-739 and Babon, *et al.*, 1999 Electrophoresis 20, 1162-1170. Other enzymes are specific for a particular functionality, such as peroxidase for a peroxide, sulfatase for a sulfate ester, pyrophosphatase for pyrophosphate, cholinesterase for acetylcholine,  $\beta$ -galactosidase for  $\beta$ -galactosidyl ethers, and the like. Thus, linkages involving peroxide groups, ester groups, anhydride groups, ether groups, as well as other functional groups may serve as the linkage for cleavage.

Instead of having a moiety resulting in complex recognition by an enzyme, one may have a cofactor bonded to the primer. By providing for cycling of the cofactor, so that the proper oxidation state of the cofactor is restored after the enzymatic reaction, the signal from a single target molecule can be amplified as described above. Various cofactors include NAD, NADP, FAD, ascorbate, glutathione, etc. and the reduced analogs thereof, when used in conjunction with metal ions such as copper or iron, as

described in Oikawa and Kawanishi, 1996 Biochemistry 35, 4584-4590; and Dreyer and Dervan, 1985 PNAS USA, 82, 968-972.

Another system employs electromagnetic radiation sensitizers, where the moiety bonded to the primer transfers energy to the linkage with resulting cleavage of the linkage. For the most part, the energy transfer moieties are dyes absorbing light in the range of about 275 to 600nm, more usually in the range of about 275 to 350nm. Of interest are intercalating dyes bound to the primer, such as porphyrins, as described in Munson and Fiel, 1992 Nucleic Acids Res. 20, 1315-1319. Illustrative dyes include texaphrin and metal complexed texaphrin (U.S. Patent nos. 5,714,328 and 5,798,491), diazapyrene (U.S. Patent no. 4,925,937), etc., where the labile bonds cleaved by energy transfer include peroxide, e.g. oxetane, alpha, beta-dicarbonyl, o-nitrobenzyloxy, hydrindanyloxy, etc. Illustrative dyes also include texaphrin and metal complexed texaphrin (U.S. Patent nos. 5,714,328 and 5,798,491), diazapyrene (U.S. Patent no. 4,925,937), etc., where the labile bonds cleaved by energy transfer include peroxide, e.g. oxetane, alpha, beta-dicarbonyl, o-nitrobenzyloxy, hydrindanyloxy, etc.

The effector molecule may be bonded to the primer nucleic acid sequence by any convenient linking group. The length of the linking group will depend on the nature of the effector molecule, the nature of the cleavable linkage, synthetic convenience, and the like. In many instances, the effector molecule may be added to the end of the synthesized nucleic acid sequence, during the synthesis of the sequence. This will usually be possible with smaller effector molecules, under about 5,000 Dal. For larger molecules, such as enzymes, the nucleic acid will usually be added to the enzyme after the nucleic acid has been synthesized. By providing for a functionality that will react with a functionality present on the enzyme. Depending on the amino acids present in the enzyme and the nature of the active site, various strategies may be used for the conjugation. Where one or more lysines are available, one may use maleic anhydride to form the maleimide, which may then be conjugated with a thiol to form the thioether. Where a cysteine is present on the enzyme, one may use active halogen, such as a benzyl halide for conjugation to form a thioether. With lysines present, one may use reductive amination with an aldehyde in the presence of a hydride supplying reductant.

Techniques for conjugating proteins to nucleic acids are well known in the literature.

See, for example, WO99/41273.

Similarly, the tag through the cleavable linker may be bonded to the probe sequence. The tags can be conveniently linked to the terminal nucleotide, either prior to 5 or subsequent to the addition of the terminal nucleotide during synthesis. Various techniques are described in the literature to provide substituents at the terminal position. Alternatively, the linkage of the synthetic sequence to the solid support can provide for a functional group that will be available upon cleavage from the solid support. The 10 functional group may then be used for linking the tag through the cleavable linkage. In this way, amides, esters, thioethers, or ethers can be employed for linking the tag and the probe sequence. Conveniently, the active functional groups will be protected during synthesis and these groups may be retained while the tag is conjugated to the sequence, after which they may be removed.

A tag molecule is one that can be differentiated from other tag molecules to allow for multiple target sequences to be analyzed in the same sample. The tag molecules are conveniently differentiated by their mobility in electrophoresis. They may have one or more bases for differentiation, where one of the bases will be mobility and a second basis may be a different detectable signal, such as electromagnetic radiation, as observed with fluorescence, chemiluminescence, electrochemiluminescence, electrical signal, etc. In this manner, the different tags will be unique as to the two or more differences, so that two tags may have the same mobility, if they may be distinguished by their fluorescence. The tags have three characteristics, which may be more or less associated with different regions of the tag. A tag will have mass, charge and be capable of detection. The mass may come from neutral molecules, charged molecules or combinations thereof. The portion of the tag associated with the mass may provide properties other than mass, such as solubility, enhancing the signal by providing for a propitious environment, allow for easy handling, synthesis or the like. Mass enhancing entities may be a single molecule without a pattern of repeating groups or have a repeating group, where the repeating group will be the individual unit of mass. The 25 repeating unit may be a divalent group, such as a hydrocarbylene or substituted hydrocarbylene, such as alkylene, e.g. methylene, ethylenene, propylene, etc., having 30

neutral substituents, such as ethers, thioethers, hydroxyl, thiols, esters, both organic and inorganic, cyano, nitro, halo, etc., may be hydrocarbylene separated by a neutral heteroatom, e.g. oxygen, sulfur, acylamido, where the nitrogen may be the intervening group or carbamyl may be the intervening group, phosphate as the triester, etc. The hydrocarbylene will be aliphatic, alicyclic, or aromatic. Alternatively, the repeating unit may be a divalent heterocyclic group, having from 1 to 3 annular heteroatoms, which are O, S, N, or P and may be combined with hydrocarbylene. The neutral mass providing entity will generally range from about 2 to 60, more usually 2 to 36 carbon atoms and from about 0 to 8 heteroatoms for the substituted hydrocarbylene or heterocyclic groups, and will usually have from 1 to 4, more usually 1 to 3 heteroatoms as a linking unit between repeating carbon containing units. The heteroatom links may be the same or different in the chain, usually the same. The groups may be polar or non-polar, preferably polar imparting water solubility.

Groups of interest may include alkylene groups, alkyleneoxy groups, alkylmethio groups, alkylenesulfone groups, glycineamide, alanineamide, serineamide, methionineamide, copolymer of ethylene diamine and succinic acid, copolymer of ethylene glycol and fumaric acid, and the like. The mass enhancing moiety will be selected based on considerations of convenience, stability under the conditions of the determination, synthetic convenience, interactions with other components of the tag and of the determination, and the like.

The charge enhancing moiety may be a single molecule or an oligomer of the same or different units, usually have the same type of linking group. Thus, the chain may be of the same nature as the mass enhancing moiety, except that the units will contribute charges to the moiety. The charges may be positive or negative. For the negative groups, illustrative groups are acyl groups, both organic and inorganic, such as carboxyl, sulfonate, sulfinate, phosphate, phosphinate, borate, phenolic hydroxyl, etc. For the positive groups, illustrative groups are amine and ammonium, sulfonium, phosphonium, metal chelates, metallocenes, etc. For the most part, oligopeptides are convenient, providing amine groups as in lysine, arginine and histidine and carboxyl groups as in aspartate and glutamate, and these will be preferred. Alternatively, alkylene imines and oligomers thereof may be used, particularly N-alkylated imines.

Depending on the nature of the moiety, there may be one or two neutral units between charged units. Generally, the number of charges will be in the range of about 1 to 10, usually 1 to 6, and frequently in the range of 1 to 4. The moiety will usually have from about 2 to 60 carbon atoms, more usually 2 to 30 carbon atoms and at least one heteroatom and may have as many as 30 heteroatoms, usually not more than about 20 heteroatoms associated with the functional groups indicated above.

The detectable moiety will generally be detectable by electromagnetic radiation, e.g. luminescent, or electrically, by its redox potential. For the former, the label may be fluorescent, chemiluminescent or phosphorescent, particularly fluorescent. For the electrically detectable label, the label will usually be an organometallic compound.

The tag will generally have a molecular weight in the range of about 100 Dal to 5 kDal, depending on the number of tags necessary for the determination. The more tags required for analysis, the larger will be the range in molecular weight.

The tags may be prepared by conventional synthetic techniques. The tag may be preprepared and attached to the nucleic acid sequence or be added to the sequence in the process of synthesizing the nucleic acid sequence on a support. Depending on the nature of the cleavable linkage, the cleavable linkage may be attached to either the nucleic acid sequence or the tag, prior to joining the tag to the nucleic acid sequence.

The target sequence may come from a variety of sources, both prokaryote and eukaryote, unicellular and multicellular organisms, mammalian and plant, etc. The target nucleic acid may be RNA or DNA, where the nucleic acid may be mRNA, tRNA, cDNA, chromosomal DNA, plastid DNA, mitochondrial DNA, etc., where the DNA may be direct from the host, reverse transcribed from RNA, amplified with PCR or other amplification scheme, or the like. Depending on the source of the DNA, the DNA may be subject to further processing, such as desalting, chromatography, separation into fractions, e.g. nucleosomes, mitochondrial DNA, nucleolar DNA, chromosomal DNA or individual chromosomes, etc. Desirably the fragments will be less than about 1cM, preferably less than about 0.1cM and usually at least about 200nt, more usually at least about 300nt, usually single stranded. The DNA will be dispersed in an appropriate hybridization medium, generally having a conventional buffer, depending upon the nature of the effector. Buffers include phosphate, borate, HEPES, MOPS, acetate, etc.

The buffer concentration will generally be in the range of about 10 to 200mM, more usually in the amount of 50 to 200mM. Other salts may be present to enhance the stringency, such as sodium chloride, sodium dodecyl sulfate, potassium chloride, magnesium chloride, etc., generally present in the range of about 0 to 200mM. For enhancing stringency, small amounts of non-interfering water-soluble organic solvents may be included, such as methanol, ethanol, acetonitrile, hexamethylphosphoramide, dimethylformamide, etc, generally being present in the aqueous solution in less than about 20 volume %, usually less than about 10 volume %. For the most part, to enhance stringency, elevated temperatures will be used, by itself or in conjunction with one or more of the other stringency-enhancing agents.

The different agents are combined with the sample in an appropriate medium as described above at a temperature, usually in the range of about 10° to 60°C, more usually in the range of about 20° to 50°C. The temperature will be selected in accordance with the melting temperatures of the primer and the probe, the stability of the enzyme, when an enzyme is present, the desired rate of cycling, the necessary residence time of the probe in order to have a reasonable rate of reaction, generally in the range of about 1min to 6h, usually in the range of about 5min to 3h and preferably in the range of about 5min to 1h. Where enzyme is labile under the conditions of the reaction, additional reagent may be added at intervals to maintain the reaction rate. Where one is only interested in a qualitative answer, it is only necessary that the reaction occur to provide a detectable amount of product and the amount of enzyme and time are not critical.

The concentration of the DNA sample will generally be in the range of about 1pM to 1mM, more usually in the range of about 1pM to 1microM. The primer and probe will generally be at least about equal to and usually in excess of the anticipated concentration range of the target DNA. The primer will be present in a mole ratio of about 1-!0:1, more usually 1-5:1, while the probe will usually be in a mole ratio of 1-500:1, usually in the range of about 10-200:1. The important factor for the primer is that most, if not all, of the target sequence present will be bound by the primer. For the probe, one will usually want an excess to allow for amplification. However, where the probe can undergo cleavage in the absence of being bound to the target DNA, the

100-006830  
15  
20  
25  
30

greater the amount of probe present, the greater the background. Therefore, the optimal amount of probe will vary with the nature of the probe and the amount of background reaction as compared to reaction when the probe is bound.

There will be at least one target sequence to be measured, more usually at least

5 2. Preferably, the number of target sequences to be measured will be at least about 10 and may be a 1,000 or more, usually not more than about 500, more usually not more than about 200. For each target sequence, one primer and one probe will be added to the assay mixture. In addition, all the reagents necessary for releasing the tag from the nucleic acid sequence will also be added. Depending on the nature of effector system,  
10 the temperature may be thermally cycled or maintained during the course of the reaction. Generally the temperature necessary for amplification of the signal will be at least about 35°C, more usually and least about 45°C and not more than about 80°C, usually not more than about 65°C. The reaction time will be sufficient for at least a sufficient amount of tag to be released for detection, usually requiring at least one minute, more usually requiring at least about 5min and not more than about 24h, usually not more than about 12h and preferably not more than about 6h.

15 After the reaction is complete, the assay mixture may be analyzed. The assay mixture may be concentrated, diluted or be unchanged when added to the capillary electrophoretic analytical device. In some instances, one may wish to separate the tags from the other material present in the assay mixture. Depending on the nature of the tags, the tags may be separated by precipitation of the tags or the nucleic acid,  
20 chromatography, extraction, or the like. The liquid medium enriched for the tags may then be used for analysis. Usually, the electrophoresis will be performed with a sieving medium that provides for clean separation of the tags. The conditions for performing the electrophoresis will vary with the nature of the tags and will be evident to those of  
25 skill in the art.

30 Alternatively, one may use liquid chromatography to separate the tags, the different mass/charge ratios providing for different mobility. The packing of the channel may provide for ion exchange as well as sieving through the column. Depending on the density of the packing various methods of providing for flow through the column may be

employed. One may use electroosmosis for a lightly packed column or pumping or applied pressure for more densely packed columns.

In order to minimize the release of the tags adventitiously from the probe, various protective mechanisms may be employed. One mechanism involves using a stem and loop structure, where the stem provides protection for the labile bond. The labile bond can be in one of the strands involved in the stem. In those situations, where cleavage is as a result of DNA-RNA hybridization, the stem may involve two RNA or two DNA strands, which would not be recognized by the nuclease. Where the cleavage is as a result of a mismatch, the two strands of the stem would be complementary.

Alternatively, where the bond is chemically scissile, the presence of the other strand would inhibit chemical attack. In the event of a photolytically labile bond, subject to cleavage in the presence of a sensitizer, the other strand would provide for light absorption without energy transfer, so as to protect the photolytically labile bond. The stem need not be very long, sufficiently long to be reasonably stable under the conditions of the assay, but resulting in a more stable binding when the loop, and potentially a portion of one strand of the stem is bound to the target sequence. By having at least a portion of the probe in the stem and loop conformation, the potential for false positives will be diminished. The particular length of the stem, its melting temperature and the choice of nucleotides may be determined empirically in light of the mechanism of cleavage, the temperature at which the determination is made, and the like.

The following examples are offered by way of illustration and not by way of limitation.

25 **EXPERIMENTAL**

**Example 1**

**Singleplex SNP Detection using a CvNY2A Nickase Reaction**

Target: . . . CGCAA TAGCC TAGCA GT(c)GC . . .

30 . . . GCGTT ATCGG ATCGT CA(g)CG . . .

Primer: 5'- CGCAATAGCCTAGCCACTTG-3'

Probe: \*ACLA001-5'-CAAGTCGTcGC -3'

^

5 \* eTag fluorescent group: Fluorescein(C<sub>6</sub>P)<sub>1</sub> (ACLA001)

^ indicates site of cleavage

c indicates base complementary to SNP locus

10 10 pmole primer is mixed with 100 ng template DNA to a final volume of 10  $\mu$ l in  
Nickase buffer composed of 10 mM Tris-HCl, pH 8.3, 10 mM MgCl<sub>2</sub>, 150 mM KCl. The  
sample is heated to 95°C for 5 min, is then annealed by cooling to 42°C for 15 min,  
followed by 15 min at 30°C. 1  $\mu$ l of annealing reaction is combined with 1 nmol of probe  
and 0.1 units of CviNY2A in Nickase buffer plus 100  $\mu$ g/ml BSA to a final volume of 10  
 $\mu$ l. The sample is incubated at 30°C for 30 min, then heated to 65°C for 10 min to  
inactivate enzyme.

15 The released electrophoretic tag is separated from undegraded probe using a  
microchannel device described in detail in PN 5,900,130, and having a configuration  
generally as sketched in FIG. 1a of PN 5,900,130. Briefly, this microchannel device is  
constructed by forming a base plate and cover of acrylic polymer, with apposing  
surfaces bonded together by thermobonding. The microchannel structure corresponds  
to two crossed linear channels of dimensions 8 mm and 5.2 cm in length. The channel  
has a trapezoidal cross-section, measuring at widest about 80  $\mu$ m and at narrowest  
about 30  $\mu$ m, with an average depth about 20  $\mu$ m. At the termini of the channels, holes  
of 2 mm in diameter are drilled as buffer reservoirs. Platinum electrodes are introduced  
20 into each of the four reservoirs.

25 The assembled device is loaded with buffer by filling reservoirs 1, 2, and 4 with  
buffer and then applying a vacuum at reservoir 3 to draw the buffer into the channels  
11, 12, 13, 14; then reservoir 3 is filled with buffer. Buffer is removed from sample  
supply reservoir 4, and replaced with 5  $\mu$ l of sample material. Voltages are applied in  
stages as indicated in the table below, yielding movement of sample as illustrated in  
30 FIGS. 1a,b and 4a,b of PN 5,900,130. The "injection" stage draws sample from  
reservoir 4 towards reservoir 2, with the borders of the sample stream constrained by

buffer drawn from reservoirs 1 and 3 towards reservoir 2. Voltages are then switched to the "separation" stage settings, causing a unit of the sample stream located in the microchannel intersection to enter channel 13, migrating towards reservoir 3, and causing remaining sample in the stream bordering the sample unit to be pulled away from the channel intersection and back towards reservoirs 2 and 4. Components of the sample are separated in channel 13 according to their electrophoretic mobilities, then detected by fluorescence using an argon-ion laser operating at 488 nm to excite the fluorescein, and a photomultiplier tube with emission filter for detection at 520 nm.

10 Running conditions:

	1	2	3	4	time	
Injection		365V	730V	0	365V	60s 100V/cm
Separation	0	300V	1420V	280V	500s	250 V/cm

15 Example 2

Multiplexed SNP Detection using an ExoIII Exonuclease Reaction

Target #1:

Target: ... CTAGCAGT(c/g) G CAGGATACTG ACCTGCGC...  
... GATCGTCA(g/c) C GTCCTATGAC TGGACGCG ...

Primer: 5'-CAAGTCGGAT ACTGACCTGC AT-3'

Probes: 5'-GT<sup>c</sup>GCCACTT G-3'-ACLA001\*

5'-GT<sup>g</sup>GCCACTT G-3'-ACLA016\*

Target #2:

Target: ... GAGATCCA(c/g)CAGCA AAGTC CGAGT CGGT ...  
... CTCTAGGT(g/c)GTCGT TTCAG GCTCA GCCA ...

Primer: 5'-CAAGTACCAA AGTCCGAGTC GAA -3'

Probes: 5'-CAcCACTACT TG-3'-ACLA018\*

^

5 5'-CAgCACTACT TG-3'-ACLA019\*

^

\* indicates eTag fluorescent groups(all containing Fluorescein):

ACLA001: FI(C<sub>6</sub>P)<sub>1</sub>-oligonucleotide

ACLA016: FI(C<sub>6</sub>P)<sub>1</sub>(C<sub>3</sub>P)<sub>1</sub>-oligonucleotide

10

ACLA018: FI(C<sub>6</sub>P)<sub>1</sub>(C<sub>3</sub>P)<sub>2</sub>-oligonucleotide

ACLA019: FI(C<sub>6</sub>P)<sub>1</sub>(C<sub>3</sub>P)<sub>3</sub>-oligonucleotide

^ indicates site of cleavage

c or g indicates base complimentary to SNP locus

15 10 pmole of each primer is mixed with 100 ng template DNA to a final volume of 10  $\mu$ l in exonuclease buffer composed of 66 mM Tris-HCl, pH 8.3 and 6.6 mM MgCl<sub>2</sub>. The sample is heated to 95°C for 5 min, then annealed by cooling to 42°C for 15 min, followed by 15 min at 30°C. 1  $\mu$ l of annealing reaction is combined with 1 nmol of each of the four probes and 0.1 units of *ExoIII* in exonuclease buffer to a final volume of 10  $\mu$ l. The sample is incubated at 30°C for 30 min, then heated to 65°C for 10 min to inactivate enzyme.

20 The released electrophoretic tag is separated from undegraded probe using a microchannel device as described in Example 1. The assembled device is loaded with buffer followed by 5  $\mu$ l of sample material as described above. Voltages are applied in stages as indicated in the table below, yielding movement of sample as illustrated in FIGS. 1a,b and 4a,b of PN 5,900,130. Components of the sample are separated in channel 13 according to their electrophoretic mobilities, then detected as described in Example 1.

25

30

Running conditions:

	1	2	3	4	time
Injection		550V	1100V	0	550V 60s
Separation	0	450V	2130V	420V	500s

It is evident from the above description and examples that the subject method provides for an efficient method for identifying nucleic acid sequences or one or a few nucleotides in a nucleic acid sequence. The methodology provides for multiplexing, so that a single sample may be interrogated as to one or numerous aspects of the sample.

5 In this way, analyses, such as single nucleotide polymorphisms, mutations, expression profiles, etc., may be analyzed in a single medium and the results determined by separating the tags. In addition, by appropriate choice of conditions, the signal may be amplified, so that a single target results in a plurality of detectable products.

All publications and patent applications mentioned in this specification are

10 indicative of the level of skill of those skilled in the art to which this invention pertains.  
All publications and patent applications set forth herein are incorporated by reference to  
the same extent as if each individual publication or patent application was specifically  
and individually indicated to be incorporate by reference.

The invention now having been fully described, it will be apparent to one of ordinary skill in the art that many changes and modifications can be made thereto without departing from the spirit or scope of the appended claims.